


BRUKER WALTHAM<sup>®</sup> the science behind SPILLERS<sup>®</sup> & WINERGY<sup>™</sup> The University of Nottingham

## PROTEOMIC IDENTIFICATION OF SECRETED BIOMARKERS IN AN EXPLANT MODEL OF EARLY OSTEOARTHRITIS

A.L. Clutterbuck, J.R. Smith, D. Allaway, P. Harris and A. Mobasher



School of Veterinary Medicine and Science

The University of Nottingham

## Introduction

- Molecular changes in early OA are poorly understood
- Comparative models are beneficial for researching OA - equine
- Proteomics;
  - increasingly being applied in OA research and in basic cartilage biology to understand the composition of the chondrocyte proteome \*
  - not extensively applied to serum free explant cultures

\* Ruiz-Ramero, C., López-Armada, M. & Blanco, F.J. (2005) *Proteomics*, 5: 3048-3059.

The University of Nottingham




## Aims and Objectives

- Can proteomics and cross-species peptide matching be applied to a serum free equine cartilage explant model to identify secreted proteins in early OA?
- The cross-species approach may help validate biomarkers identified in studies using experimental animals and human cartilage samples

School of Veterinary Medicine and Science

The University of Nottingham

## Method- Explant Model


- Distal forelimbs of 3 horses euthanized for reasons other than research
 
  - Articular cartilage obtained from MCP joints and washed
 
- Full depth cartilage discs excised and placed into the wells of a 24 well plate containing serum free DMEM + 2% antibiotic
 

School of Veterinary Medicine and Science

The University of Nottingham

## Method- Study Design

- After overnight incubation at 37°C, the media was removed and replaced
- Explants were incubated at 37°C for 5 days with;



School of Veterinary Medicine and Science

The University of Nottingham

## Sample Preparation


- Reduction of soluble proteins - DTT
- Blocking the thiol groups - iodoacetamide
- Protein precipitation - ice-cold acetone
- Protein quantification - Lowry Assay
- Protein resuspension - ammonium bicarbonate (at 2µg/µl)
- Peptide fragmentation - trypsin
- Termination of trypsin digestion - formic acid

School of Veterinary Medicine and Science

## Separation & Detection




- Separation
  - on a C18 pepmap column using a Dionex U3000 chromatography platform
- Detection
  - Bruker HCT PTM discovery ion trap instrument.




School of Veterinary Medicine and Science

## MS/MS Ion Search



School of Veterinary Medicine and Science

## Chromatographic Profiles




School of Veterinary Medicine and Science

## MS Scan



- The 4 most abundant ions in each MS scan were selected for fragmentation.




School of Veterinary Medicine and Science

## Sequence Matching



- Fragment patterns for each peptide were compared to the mammalian entries in the SwissProt database using MASCOT search engine



Mascot Search <http://www.matrixscience.com/>

## Results- Cross Species Matching



Protein	Bovine	Dog	Chimp	Horse	Human	Mouse	Pig	Rabbit	Rat	Sheep
FMOD	*					*			*	
CILP-1					*	*	*			
CHAD	*				*	*			*	
MMP-1					*	*				
PRELP					*	*				
PGS-1	*			*						
PGS-2	*	*		*	*			*	*	
TR11B					*	*				

## Results- Consistently Identified Proteins




- The most abundant proteins secreted into the media were identified
  - Proteins constitutively expressed by cartilage
  - Found in other proteomic studies \*

Protein	Control (18)	IL-1 $\beta$ (18)
Aggrecan Core Protein	18	18
Fibronectin	17	18
Cartilage Oligomeric Matrix Protein	18	17
Clusterin	18	16

- Validate the model

\* Stevens, A.L., Wishnok, J.S., Chai, D.H., Grodzinsky, A.J. & Tannenbaum, S.R. (2008) *Arthritis & Rheumatism*, 58 (2): 489-500

## Differentially secreted proteins- IL-1 $\beta$ specific



School of Veterinary Medicine and Science

## MMP-3 -horse



- Increased activity in SF of equine OA joints\*
- Increased in equine chondrocyte cultures stimulated with IL-1 $\beta$  \*\*




\* Brama, P.A. et al. (2000) *Annals of the Rheumatic Diseases*, 59(2): 155-157.

\*\* Tung, J.T. et al. (2002) *Canadian Journal of Veterinary Research*, 66(1):19-25.

## TSP-1 -bovine, human, mouse



- Adhesive glycoprotein that mediates cell-to-cell and cell-to-matrix interactions.
- Increases in early human OA cartilage lesions\*



\* Pfander, D., et al. (2000) *Annals of the Rheumatic Diseases*, 59: 448-454.

## Conclusions



- First proteomic study of a serum free equine cartilage explant model
- Identified several proteins with well-established ECM functions
- Other proteomic studies have identified some of the proteins being up-regulated in early OA thus validating the model
- Identified proteins that could be potential biomarkers in the supernatants of IL-1 $\beta$ -stimulated equine cartilage explants
- Some of the proteins identified may participate in early repair responses and may be useful for the detection of early changes in OA

School of Veterinary Medicine and Science

## Acknowledgements



This work was supported by the Biotechnology and Biological Sciences Research Council (BBSRC, grant BS/S/M/2006/13141) and the WALTHAM Centre for Pet Nutrition



School of Veterinary Medicine and Science

Thank you

